

Passwordless troubleshooting_Updated

Article ID: KB0024121

Created: 2023-08-22 19:42:37

Updated: 2023-08-22 19:43:07

Author: {'link':

'https://fmcnadev.service-now.com/api/now/table/sys_user/31a85ef41bb989501d67db1dcd4bcb92',

'value': '31a85ef41bb989501d67db1dcd4bcb92'}

Category: {'link':

'https://fmcnadev.service-now.com/api/now/table/kb_category/c2c873fd1b7cb91070dacao3604bcb79',

'value': 'c2c873fd1b7cb91070dacao3604bcb79'}

Passwordless Project

Troubleshooting Steps

Problem: Scanner lights don't blink when a Yubikey is scanned, and scanning is not working.

Solution: Reseat the scanner USB cable. If the scanner still doesn't work, plug in a mouse to the USB port and verify the mouse works. Inform the Biomed if there is a problem with the USB.

Problem: Chairside shows a red scanner icon and scanning is not possible.

Solution: Restart the T machine. If the scanner icon remains red inform the Passwordless team. For Biomed open a terminal and with the fresenius user execute the command: curl localhost:8282 .You should receive a confirmation message the scanner is working.

```
{"service": "Passwordless Python Daemon", "version": "1.0.0-August-20210817-V1",  
"smart_card_reader_status": "Found"}
```

Problem: After registering a Yubikey the user gets an error message when scanning.

Solution:

Step 1: Please verify User Account Status in OKTA. It should be ACTIVE & YubiKey should be registered in OKTA . Use OKTA Admin login to verify the YubiKey Status. It should be ACTIVE.

Step 2: After checking step1, Verify whether user is able to login with the below given URL's

OKTA

<https://fmc.okta.com/>

OR

Service Now

<https://fmcna.service-now.com/>

Step 3: Above given step1 & 2 are fine, still user is having the below given issue, please reach out to Passwordless Support Group Chat and share the Emp id & email id , Serial number of the Yubikey. Keep the INC in your L1 Queue. Based upon IAM /Chairside Team resolution, Help Desk Team can call back and ask users to try the resolution.

Problem: Invalid Username and/or Password error in chairside when try with Manual option.

Step 1: Please verify User Account Status in OKTA. It should be ACTIVE. If it's LOCKED, Help Desk team can UNLOCK by using OKTA Admin tool.

