Adding User to Tableau Server

Article ID: KB0039498

Created: 2024-12-19 14:01:54 **Updated:** 2024-12-19 14:01:54

Author: {'link':

'https://fmcnadev.service-now.com/api/now/table/sys_user/acaf843397d74958f7e3bb8fe153afe3',

'value': 'acaf843397d74958f7e3bb8fe153afe3'}

Category: {'link':

'https://fmcnadev.service-now.com/api/now/table/kb_category/2066c40cc3821a100acd3300150131d1',

'value': '2066c40cc3821a100acd3300150131d1'}

1. Login to Tableau Server Qual: https://qdata-qual.intra.fmeone.com Prod: https://qdata.intra.fmeone.com

2. Select the site to which the user needs access.

Click on Add User

Click on Active Directory User

Provide domain name\username

Select the required Site role

Click on Import Users